


DX75480 LAPPING MACHINE

***(INDISPENSABLE MICROSCOPE DX75680)**

**** (RECOMMENDED MICROSCOPE LIGHT DX75681)**





This lapping machine has to be used for reconditioning the ball seat internal valve cone on the CR. Bosch type injectors, in order to re-establish the ball tightness in its seat.





MAIN CHARACTERISTICS:

Electric power supply: 230 V - 50 Hz
 ON/OFF: It goes On pushing and Off releasing the punch.

KIT CONTENTS

<p>A </p> <p>B </p> <p>C </p>	<p>Three differents spindles with different angled cones: A: DX 75487 B: DX 75485 C: DX 75486 (technical informations are available on our website codec.ditex.it)</p>
	<p>DX75483 Rotating valve holder to insert the valve for CRIP1 - CR12.1 - CR12.2 - CRIN - CRIN2</p>

OPTIONAL






	<p>DX 75488 Universal spindle with more acute angle to re-establish the ball seat in depth</p>
	<p>DX75482 Rotating valve holder to insert the valve for CRIN3</p>

ACCESSORIES (NOT INCLUDED IN THE KIT)

GRINDING STRIPS FOR LAPPING MACHINE

Pack of 50 strips in pre-cut abrasive paper in 5 different graduations



-  DX75490 (raw)
-  DX75491 (medium)
-  DX75492 (fine)
-  DX75493 (polishing)
-  DX75494 (very raw)

DX75470 KIT TO GRIND SETTING SHIMS




***(INDISPENSABLE DIGITAL MICROMETER GAUGE DX75391 WITH HOLDER DX75392)**
**** (RECOMMENDED SYNTHETIC LIQUID DX75394)**



PRE-CUT ABRASIVE PAPER FOR DX75470

Pack of 10 pre-cut abrasive paper in 3 different graduations



-  DX75723 (fine 500)
-  DX75724 (medium 320)
-  DX75729 (raw 180)